

Contec® NeutraKlean

Neutral Detergent with Purified Water

Contec NeutraKlean is a gentle and non-foaming detergent designed for the cleaning of life science cleanrooms.

With a pH after irradiation between 6.5 and 7.5, Contec NeutraKlean is ideal for carrying out the low-level cleaning required after maintenance, and to remove disinfectant residues and general soil. The use of a surfactant/detergent helps to break the binding of particles and dirt to a surface and allows them to be more easily removed.

The product is 0.2 µm filtered and filled under ISO Class 5 airflow and packaged in an ISO Class 7 cleanroom. Ready to use, the product is available in a 1L trigger spray with a “bag in bottle” system to protect the contents during use. A non-sterile version is available in a 5L capped container.

Both the sterile and filtered product are provided double bagged to aid the transfer of the product into a pharmaceutical cleanroom. The linear tear bags are easy to open even when wearing gloves.



Filled under ISO Class 5 air and packaged in an ISO Class 7 cleanroom, Contec NeutraKlean is available in a 1L trigger spray with a “bag in bottle” system to protect the contents during use.

FEATURES	BENEFITS
Neutral pH	• Suitable for use in all areas of cleanroom as gentle on surfaces
Filtered to 0.2 microns under ISO Class 5 airflow in a ISO Class 7 cleanroom	• Ensures the product is free from contamination and particulates
Low foaming detergent action	• Ideal for use on floors, walls, and benches
Low level of surfactant	• Low residue which can be easily removed
Trigger spray can be set to jet or spray	• Large droplet size reduces the risk of inhalation and provides good surface coverage
Trigger spray and “bag in bottle” protected system	• Prevents bottle contents from becoming contaminated during use
Double bagged packed in linear tear packaging	• Each bag is easy to open even when wearing gloves • Facilitates transfer of cleaning products into isolators and RABS

Full documentation available upon request.

Part No.	Description	Quantity
SBT103NK	Contec Irradiated NeutraKlean, 34 oz. (1L), w/trigger spray, sterile	6 bottles/case
FBC503NK	Contec NeutraKlean, 1.3 gallon (5L)	2 bottles/case

Visit our website for more information: www.contecinc.com.

Copyright © 2017 Contec, Inc. All rights reserved. DIS032 062017

Contec, Inc
Spartanburg, SC
United States
Tel: +1 (864) 503-8333

Contec
Vannes Cedex
France
Tel: +33(0) 2 97 43 76 98

Contec Cleanroom Technology (Suzhou) Co. Ltd
Suzhou
China
Tel: +86-512-6274 4050

